


# Peste des Petits Ruminants: Diagnosis training

Olivier KWIATEK

CIRAD

Research Unit 15: Control of emergent and exotic animal  
diseases

Montpellier - France



**cirad**

LA RECHERCHE AGRONOMIQUE  
POUR LE DÉVELOPPEMENT

# Peste des Petits Ruminants diagnosis

## Part 1

## ELISA

# ELISA

→ **ELISA = Enzyme Linked ImmunoSorbent Assay**

- **The Enzyme Linked ImmunoSorbent Assay is an immunoenzymatic test of detection**
- **The reaction antigen/antibody is visualized by a colored reaction of a substrate with an enzyme linked to the antibody**


# Materials and equipment required



**Spectrophotometer**



**Multichannel pipette**



**96- well polystyrene plate**



**micropipettes**

# Important points for ELISA test

- **Pipetting** → the pipette must be regularly verified
- **Washing** → use high quality purified water for preparation of reagents
  - be careful for the cleaning of the material
- **Temperature** → incubate at the correct temperature ( $37 \pm 5$ )°C or ( $21 \pm 5$ )°C
- **Timing** → use a timer
- **Spectrophotometer** → use the correct wavelength

# PPR c-ELISA

**IDVET**  
innovative  
diagnostics

## ID Screen<sup>®</sup> PPR Competition


Kit for the detection of antibodies against the Peste des Petits Ruminants  
(PPR) virus in sheep and goat serum and plasma by competitive screening  
ELISA

For *In vitro* use



PPRC ver 1209 GB

# PPR c-ELISA



# Testing procedure

- ❖ Allow all the reagents to come to room temperature ( $21 \pm 5$ )°C before use
- ❖ Homogenize all reagents by inversion or vortex

1. add 25  $\mu$ L of dilution buffer 13 to each well
2. add 25  $\mu$ L of the positive control to well A1 and B1
3. add 25  $\mu$ L of negative control to well C1 and D1
4. add 25  $\mu$ L of each sample to be tested to the remaining wells
5. Incubate 45 min at 37°C
6. Wash each well 3 times with approximately 300  $\mu$ L of Wash solution. Avoid drying of the wells between washings
7. Prepare the conjugate 1X by diluting the conjugate 10X to 1/10 in dilution buffer 4
8. Add 100  $\mu$ L of the conjugate 1X to each well
9. Incubate 30 min at 21°C
10. Wash each well 3 times with approximately 300  $\mu$ L of Wash solution. Avoid drying of the wells between washings
11. Add 100  $\mu$ L of the substrate solution to each well
12. Incubate 15 min at 21°C in the dark
13. Add 100  $\mu$ L of stop solution to each well in order to stop the reaction
14. Read and record the O.D. at 450 nm

# Validation of the test

## The test is validated if:

- The mean value of the negative control O.D. is greater than 0.7

$$\rightarrow OD_{NC} > 0.700$$

- The mean value of the positive control O.D. is less than 30% of the  $OD_{NC}$

$$\rightarrow OD_{PC} / OD_{NC} < 0.3$$

# Interpretation

- For each sample, calculate the competition percentage:

$$S/N = (OD_{\text{sample}}/OD_{\text{NC}}) \times 100$$



Result	Status
$S/N \leq 50\%$	<b>POSITIVE</b>
$50\% < S/N \leq 60\%$	<b>DOUBTFUL</b>
$S/N > 60\%$	<b>NEGATIVE</b>

# Peste des Petits Ruminants diagnosis

## Part 2

### PCR

# PPR diagnosis by RT – PCR (1)

→ **RT = Reverse Transcription**

- to synthesize cDNA from RNA samples

→ **PCR = Polymerase Chain Reaction**

- to amplify DNA


Pcr.exe

→ **RT - PCR = Reverse Transcription + Polymerase Chain Reaction**

- to synthesize cDNA from RNA samples and then amplify DNA


# PCR : Polymerase Chain Reaction

30 - 40 cycles of 3 steps :



Step 1 : denaturation


1 minut 94 °C



Step 2 : annealing

45 seconds 54 °C

forward and reverse primers !!!



Step 3 : extension

2 minutes 72 °C

only dNTP's

Reaction tube content:

-DNA to amplify : 50ng to 500ng

-2 sépcifics primers of the targed DNA

-dNTP

-Taq Polymerase

Final volume : 20  $\mu$ l to 100  $\mu$ l

Equipment : thermocycler

# PPR diagnosis by RT – PCR (2)

→ 3 main steps

1. extraction

2. amplification

3. analysis




RNA extraction from samples

cDNA synthesis (RT)




One-step RT - PCR



Electrophoresis on 1.5% agarose gel

# Limits of this technique

- need results visualization on agarose gel
- need use of Ethium Bromide (= risk)
- could be a source of contamination for the laboratory
- difficulties to record data
- qualitative results only (- or +)
- need 4 to 8 hours to obtain final results for tested samples




**Development of Real-Time PCR  
protocols**

# PPR one-step RT-PCR protocol

- Extract RNA from samples using 150 µl of the viral suspension
- One-step RT-PCR

		One step Qiagen Mix 1 Reac (50µl)	Mix10 reactio ns	Mix 20 reactio ns
RNA		5	5	5
5X Tp qiagen		10	100	200
dNTPmix qiagen		2	20	40
Q solution		10	100	200
NP3 10uM		3	30	60
NP4 10uM		3	30	60
Qiagen enz mix		2	20	40
H2O		15	150	300
Programme				
RT		30 min 50°C		
Denaturation		15 min 95°C		
40 cycles		30 sec 94°C		
		30 sec 60°C		
		1 min 72°C		

# PCR versus Real Time PCR



---


# The Realtime PCR

---

*« Quantitative PCR », « Real time PCR », « Kinetic PCR »,  
« Real time RT PCR » (sur ADNc), « FRET PCR » (Förster resonance  
energy transfer )*

- PCR product formation followed over the time by in situ fluorescence emission (many possible approaches).
- Determination of the number of cycles from which the PCR product is detectable : threshold cycle number = "Cycle threshold" (or "crossing point" = Ct.
- Time of appearance of the threshold signal (= Ct) depending on the quantity of matrix initially present in the amplified sample.  
Linearity over 5 to 7 orders of magnitude.

# Monitoring de l'amplification



---

Used for quantify:

- **Genomic DNA** : for example, detection/quantification of microorganisms in tissue
- Or **cDNA** produce from mRNA after reverse transcription

Quantitative PCR relies on continuous measurement of fluorescence release during PCR cycles reflecting PCR product accumulation.

---

## Main techniques of fluorescence emission during PCR

We distinguish several approaches of QPCR according to the mechanisms leading to the increase of the fluorescence during the amplification:



- Fluorescence emission by an intercalating agent (e.g. SybrGreen® I) inserted into double-stranded DNA
- Fluorescence release based on :
  - The 5'-3' exonuclease activity of the Taq polymerase exerted on a probe (eg TaqMan®), the degradation of which during the elongation of the PCR fragment makes it possible to remove the inhibitory action of a quencher on a fluorophore

# Génération de la fluorescence (1)


Avant amplification

Après amplification

SYBR Green



# TAQMAN PROBES <sup>®</sup>



## Fluorochromes *reporters* nature:

- FAM : blue
  - VIC : green
  - NED : black
  - PET : red
- Best quenching

## Advantages of Q-PCR with probes:

- Very specific (this PCR requires the use of two primers and one or more specific probes)
- Possibility of multiplexing i.e. several PCRs in the same tube: detection of different fluorochromes

**Limit :** Cost of probes (especially if many genes to quantify)

## Two technics of quantification :

-The absolute quantification: it consists in the determination of the exact number of copies of target matrix. This approach requires a calibrated range of which the exact number of copies (e.g., a range realized using a perfectly defined quantity of a plasmid carrying the gene of interest) is known.

-The relative quantification: it consists in the direct comparison of two given conditions.

Example: in the determination of the ratio of the target quantity between the two conditions (in this case the ratio of the amount of *amiF* transcripts after culture at pH 7 and at pH 5 , 5).

This approach requires the quantification in parallel for each condition of the amount of a reference gene. This reference is sometimes tricky to choose. In fact, this must be a gene whose expression does not vary between the two conditions compared.

The genes most often used as references are so-called "housekeeping" genes such as bacterial *RNA16S*, *gapdh* (encoding glyceraldehyde-3-phosphate dehydrogenase) or the gene encoding eukaryotic actin


## Standard curve technique

- Draw the standard curve :  
 $y = Ct$   
 $x = \log(\text{concentration})$
- Calculate a and b such as :  $y = ax + b$
- Report on the curve the ct values of your samples

# $C_T$ value

- $C_T$  “cycle treshold”
- $C_T$  is inversely proportional to the initial copy number
- The  $C_T$  value can be directly correlated to the starting target concentration.

RT-PCR Amplification Plot



RT-PCR Standard Curve

